EN Instructions for Use ZH 面罩
























**JOYCEeasy**  
**JOYCEeasy Full Face**  
**JOYCEeasy Full Face NV**

Nasal Mask  
Full Face Mask

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**LÖWENSTEIN**  
medical

**English**

**4**

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16

# Table of Contents

<b>1 Introduction</b> .....	<b>5</b>	<b>9 Technical data</b> .....	<b>13</b>
1.1 Model .....	5	9.1 Materials .....	13
1.2 Intended use .....	5	9.2 Characteristic pressure/ flow curve .....	14
1.3 Contraindications .....	5	<b>10 Markings and symbols</b> ..	<b>14</b>
1.4 Side effects .....	5	<b>11 Warranty</b> .....	<b>15</b>
1.5 Structure and components .....	5	<b>12 Declaration of conformity</b>	<b>15</b>
<b>2 Safety</b> .....	<b>5</b>		
2.1 Safety information .....	6		
2.2 General information .....	7		
<b>3 Product description</b> .....	<b>7</b>		
3.1 Overview .....	7		
3.2 Compatible devices .....	7		
3.3 Exhalation system .....	8		
3.4 Anti-asphyxia valve .....	8		
3.5 Pressure measurement port .....	8		
<b>4 Operation</b> .....	<b>8</b>		
4.1 Applying the mask .....	8		
4.2 Adjusting the mask .....	8		
4.3 Removing the mask .....	9		
4.4 Dismantling the mask .....	9		
4.5 Assembling the mask .....	9		
<b>5 Hygiene preparation</b> .....	<b>10</b>		
5.1 Cleaning the mask .....	10		
5.2 Tips .....	11		
5.3 Hygiene treatment (clinical sphere) .....	12		
<b>6 Disposal</b> .....	<b>12</b>		
<b>7 Faults</b> .....	<b>12</b>		
<b>8 Replacement parts</b> .....	<b>12</b>		

# 1 Introduction

## 1.1 Model

JOYCEeasy, JOYCEeasy Full Face

## 1.2 Intended use

JOYCEeasy nasal masks and JOYCEeasy Full Face masks are used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

## 1.3 Contraindications

The mask may not be used in the following situations: Immediate intubation required; loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: Pressure points and acute injuries to the skin of the face; skin allergies involving the face; deformities of the face or nasopharynx; acute pain affecting the face; cough reflex restricted or absent; claustrophobia; acute nausea.

If you are not certain whether one of these situations applies to you, please consult your attending physician or medical advisor. Please also observe the contraindications in the instructions for use of your therapy device.

## 1.4 Side effects

The following side effects may occur with use of the mask: nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

Should such side effects occur, please contact your attending physician or medical advisor.

## 1.5 Structure and components

The mask consists of a headgear, forehead cushion, forehead support, lock, retaining ring, elbow, plug for pressure measurement port, pressure measurement port, rotating sleeve, mask body, headgear clip, mask cushion and an anti-asphyxia valve (on full face masks only). The product is non-sterile and reusable.

# 2 Safety

Please read these instructions carefully. They are a constituent part of the mask and must be available at all times. In accordance with MDR (EU) 2017/745, observe the points below.

## 2.1 Safety information

### 2.1.1 Safety information

#### **Risk of injury from mask parts breaking off!**

Deteriorated mask parts or those under severe strain may come off and put the patient at risk.

- ⇒ Note useful life.
- ⇒ Check mask parts regularly and replace prematurely if necessary.

#### **Risk of injury from excessive leaking!**

Excessive leaking may lead to under-supply to the patient.

- ⇒ Activate low pressure/leak alarms on the therapy device.
- ⇒ Use the correct mask size and check that it is securely in position.

#### **Risk of injury from the supply of oxygen!**

Oxygen can deposit on clothing, bed linen and hair. In conjunction with smoking, naked flames or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Do not use electrically conductive tubes (exception: breathing tubes which can be heated electrically).

#### **Risk of injury from oxygen therapy!**

Oxygen therapies can lead to side effects in case of incorrect use or dosage.

- ⇒ Only administer oxygen as prescribed by a physician.
- ⇒ Follow the instructions for use of the oxygen system and therapy device.

#### **Risk of injury from re-inhaling CO<sub>2</sub>!**

If the mask is used incorrectly, CO<sub>2</sub> may be re-inhaled.

- ⇒ Do not close the openings of the mask.

⇒ Only apply the mask for a prolonged period if the therapy device is running.

⇒ Only use the mask within the quoted therapy pressure range.

⇒ Do not use the mask on patients who cannot take it off themselves.

⇒ Monitor patients with obstructive and restrictive lung diseases individually when using JOYCEasy masks.

#### **Risk of injury from slipping of the mask!**

If the mask slips or falls off, the therapy is ineffective.

- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.
- ⇒ Monitor patients with obstructive and restrictive lung diseases individually.

#### **Patient endangerment through application of mask during anesthesia!**

When the nasal mask is used during anesthesia, the anesthetic gas supplied can escape through the mouth and endanger the patient.

- ⇒ Do not use nasal mask during anesthesia.

#### **Risk of injury from lack of cleaning!**

The mask may show contamination which may put the patient at risk.

- ⇒ Clean the mask before using for the first time (see section entitled "Cleaning and hygiene treatment").
- ⇒ Clean mask regularly.

## 2.2 General information

- Use of third-party products may lead to functional failures and restricted fitness for purpose. Biocompatibility may also be compromised. Please note that in these cases, any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions are used.
- To prevent infection, bacterial contamination or functional impairments, follow the instructions in the section "Hygiene preparation" (see "5 Hygiene preparation", page 10).
- In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

## 3 Product description

### 3.1 Overview

Illustrations of the individual components can be found on the enclosed folding sheet (figure A1).

#### Key

- 1 Headgear
- 2 Forehead cushion
- 3 Forehead support
- 4 Lock
- 5 Retaining ring
- 6 Elbow

- 7 Plug for pressure measurement port
- 8 Pressure measurement port
- 9 Rotating sleeve
- 10 Mask body
- 11 Headgear clip
- 12 Mask cushion
- 13 Anti-asphyxia valve (on full-face masks only, see figure F)

### 3.2 Compatible devices

The mask can be used with all therapy devices which are not for life-sustaining treatment of patients with sleep apnea or ventilatory insufficiency.

With many combinations of device, the actual pressure in the mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to therapy pressure. The therapy pressure required may vary between different mask types, so prescription of a suitable therapy pressure should in each case involve adjusting/adapting therapy to suit the mask type which is also going to be used during therapy itself.

#### Quick-release cord (optional)

You can use a quick-release cord (accessory) to open the headgear in emergency situations (e.g. vomiting).

The quick-release cord consists of a special clip with a cord which you can fit to the mask instead of one of the bottom headgear clips.

In an emergency, simply pull the quick-release cord. The special clip will detach itself from the mask. The headgear is now open and the mask can be removed (**figure G**).

### 3.3 Exhalation system

The mask has an integrated exhalation system. The retaining ring and mask body are shaped in such a way that a gap forms between these parts. The exhaled air can escape through this gap (**figure A2**).

### 3.4 Anti-asphyxia valve

#### **WARNING**

##### **Risk of asphyxia if anti-asphyxia valve not working properly!**

Residues may block the mask and jeopardize the success of therapy.

⇒ Check before each use that the anti-asphyxia valve is working properly.

⇒ Clean contact surfaces between the elbow and the rotating sleeve thoroughly.

##### **Safety function of the anti-asphyxia valve**

The anti-asphyxia valve has an opening to the atmosphere through which the patient can inhale ambient air if the therapy device fails. This significantly reduces the risk of breathing CO<sub>2</sub> back in and thus also the risk of asphyxia. Inside the anti-asphyxia valve is a membrane which can adopt one of two positions:

- as long as the therapy device is in operation, the air flow forces the membrane against the opening to the atmosphere and seals it completely. The patient breathes through the therapy device.
- if the therapy device fails or if the device is not connected, the membrane is in the rest position, in other words the opening to the atmosphere is free. The patient breathes ambient air through this opening.

### 3.5 Pressure measurement port

The mask has a connection **8** which can be used for measuring pressure or supplying oxygen. If you are not using the pressure measurement port, close it with the plug **7** to reach therapy pressure.

## 4 Operation

### 4.1 Applying the mask

1. Pull the headgear **1** over the head and place the mask on the face (**figure B1**).
2. Insert headgear clips **11** (**figure B2**).

### 4.2 Adjusting the mask

1. Undo lock **4** (**figure C1**).
2. Pull the forehead support **3** in or out of the lock **4** bit by bit as necessary (**figure C1**).

Remember the most convenient position. To do this, use the scale between the forehead support and the lock (figure C1).

3. Applying the mask (see "4.1 Applying the mask", page 8).
4. Adjust headgear **1** with the aid of the hook-and-loop fasteners so that the mask seals on the face but is not too tight (figure C2).
5. To adjust the hook-and-loop fastener on top of the head undo the hook-and-loop fastener, pull on the ends of the headgear with both hands and close the fastener again (figure C3).
6. Turn retaining ring **5** counterclockwise and remove it (figure D6).
7. Remove elbow **6** from retaining ring **5** (figure D7).
8. Take rotating sleeve **9** off elbow **6** (figure D7).
9. For full-face masks only: carefully take the anti-asphyxia valve off the elbow. In the process, please note: do not damage the membrane on the inside of the anti-asphyxia valve.

### 4.3 Removing the mask

1. Undo **11** the headgear clips.
2. Slip the headgear **1** over the head and remove the mask.

### 4.4 Dismantling the mask

1. Undo the headgear clips **11** (figure D1).
2. Release hook-and-loop fasteners from headgear clips **11** and forehead support **3** and remove the headgear **1** (figure D1 + D2).
3. Remove the plug for the pressure measurement port **7** (figure D3).
4. Remove forehead cushion **2** from forehead support **3** (figure D4).
5. Take mask cushion **12** off mask body **10** (figure D5).

### 4.5 Assembling the mask

1. Push rotating sleeve **9** onto elbow **6** (figure E1).
2. Insert elbow **6** in the retaining ring **5** (figure E1).
3. Insert retaining ring in the mask body with position **5** at 11 o'clock and turn clockwise to a 12 o'clock position until it locks into place (figure E2).
4. Fit mask cushion **12** onto mask body **10** (figure E3).
5. Pull forehead cushion **2** over the forehead support **3** (figure E4).
6. Fit the plug for the pressure measurement port **7** (figure E5).
7. Attach the headgear **1** with the headgear clips using the lower hook-and-loop fasteners **11** (figure E6).
8. Attach the headgear **1** to the mask (figure E7).

9. For full-face masks only: carefully push anti-asphyxia valve onto the elbow until it audibly engages on both sides of the elbow (figure F).

In the process, please note: do not damage or trap the membrane on the inside of the anti-asphyxia valve.

You can tell whether the headgear is correctly attached by the hook-and-loop tabs. The hook-and-loop tabs on the forehead have the same color as the inside of the head bands.

## 5 Hygiene preparation

### WARNING

#### **Risk of injury from inadequate cleaning!**

Residues may congest the mask, impair the integrated exhalation system and jeopardize therapy success.

⇒ For patients with a compromised immune system or particular background of illness, disinfect mask components daily following consultation with the physician.

⇒ Thoroughly clean the contact surfaces between the elbow and the rotating sleeve.

- **This product may contain disposable items. Disposable items are intended for single use.** Therefore, use these items only once and do **not** reprocess them. Reprocessing disposable items may impair functionality and safety of the product and lead to unforeseeable reactions due to ageing, embrittlement, wear, thermal stress, the effects of chemical processes etc.

## 5.1 Cleaning the mask

1. Wash your hands before cleaning.
2. Dismantle the mask (see "4.4 Dismantling the mask", page 9).
3. Clean mask in accordance with the table below:

Action	Daily	Weekly
Wash mask parts by hand in hot water (approx. 30 °C) and with mild detergent (1 ml detergent to 1 l water) for 15 minutes.	X	
Clean mask parts (exception: mask cushion, anti-asphyxia valve and forehead cushion), thoroughly with a cloth or a soft brush for at least 3 minutes when washing. Or: Put mask parts in the top basket of the dishwasher. Select a washing program with a max. temperature of 70 °C and a max. time of 90 minutes. Use a mild dishwasher detergent (1 ml dishwasher detergent to 1 l water). Only wash mask parts in a dishwasher separately, do not mix with dirty crockery.		X
Wash the headgear by hand in hot water (approx. 30 °C) and with mild detergent (1 ml detergent to 1 l water) for 15 minutes.		X

4. Rinse all parts with clear water.

### NOTICE

#### Material damage due to abrasion!

Rubbing the mask cushion and forehead cushion can damage the coating.  
⇒ Do not rub the mask cushion and forehead cushion dry.

### NOTICE

#### Material damage due to careless use!

Careless handling can damage the anti-asphyxia valve.  
⇒ Always handle the anti-asphyxia valve with care.

5. Allow all parts to air-dry.
6. Perform a visual inspection.
7. Replace damaged parts if necessary.
8. Assembling the mask (see "4.5 Assembling the mask", page 9).

Discolorations of mask components do not impair the functionality of the mask.

## 5.2 Tips

- You can wash all parts (exception: anti-asphyxia valve and headgear) in a dishwasher at 65 °C.
- Do not dry parts in direct sunlight.
- Wash headgear before first use, as the color may run.
- Do not iron the headgear in order to keep the hook-and-loop fasteners intact.
- Do not dry the headgear in a tumble-dryer.

## 5.3 Hygiene treatment (clinical sphere)

On change of patient, inadequate hygiene treatment may lead to a risk of infection for the patient. In the event of a change of patient, subject the mask to a hygiene treatment in line with the “Information on hygiene treatment” brochure. The brochure can be found on the manufacturer’s website. We will send you this brochure on request.

## 6 Disposal

You can dispose of all parts in domestic waste.

## 7 Faults

Fault	Cause	Remedy
Facial pressure pain.	Mask fits too tightly.	Loosen headgear. Adjust position of forehead support.
Draft in the eyes.	Mask fits too loosely.	Tighten headgear. Adjust position of forehead support.
	Mask does not fit.	Use different mask size. Contact your specialist retailer.

Fault	Cause	Remedy
Therapy pressure is not reached.	Mask is not adjusted correctly.	Re-adjust mask (see 4.2, p. 8).
	Mask cushion is damaged.	Replace mask cushion.
	Tube system is leaking.	Check connection and fit of tubes.
	Air is escaping from pressure measurement connection.	Close pressure measurement port with plug (included in scope of supply).
	Anti-asphyxia valve incorrectly fitted.	Fit anti-asphyxia valve correctly (see 4.5, p. 9).
	Anti-asphyxia valve defective	Replace anti-asphyxia valve

## 8 Replacement parts

A current list of replacement parts can be ordered on the internet site of the manufacturer or through your authorized specialist dealer.

## 9 Technical data

	Nasal mask	Full-face mask
Product class according to directive MDR (EU) 2017/745	IIa	
Dimensions (W x H x D)	approx. 80 mm x 130 mm x 90 mm <sup>1</sup>	approx. 120 mm x 150 mm x 110 mm <sup>1</sup>
Weight	approx. 72 g <sup>1</sup>	approx. 103 g <sup>1</sup>
Therapy pressure range	4 hPa - 30 hPa	
Hose connection: Cone according to EN ISO 5356-1	Ø 22 mm (male)	
Pressure measurement port	Ø 4 mm	
Wide headgear clips	max. 20 mm	
Temperature range: Operation Storage	+5 °C to + 40 °C -20 °C to +70 °C	
Flow resistance at 50 l/min	0.16 hPa	0.13 hPa
Flow resistance at 100 l/min	0.64 hPa	0.53 hPa
Anti-asphyxia valve flow resistance		
• Inspiration at 50 l/min:	-	0.5 hPa
• Exhalation at 50 l/min:		0.5 hPa

	Nasal mask	Full-face mask
Switching pressure for anti-asphyxia valve	-	
• Open:		≤ 1 hPa
• Close:		≤ 2 hPa
Quoted two-figure noise emission value according to ISO 4871:		
-Sound pressure level	18 dB(A)	
-Sound power level, uncertainty	23 dB(A)	
factor:	3 dB(A)	
Service life	5 years <sup>2</sup>	
Standards applied	EN ISO 17510:2020	

<sup>1</sup> depending on size and variant.

<sup>2</sup> Useful life up to 12 months. The materials of the mask deteriorate if exposed to e.g. aggressive detergents. In individual cases, therefore, it may be necessary to replace mask parts sooner.

# CE 0197

## 9.1 Materials


Mask part	Material
Mask cushion, forehead cushion	Silicone, silicone coating
Plug for pressure measurement port	Silicone
Retaining ring (hard component)	PP (polypropylene)

Mask part	Material
Retaining ring (soft component)	TPE (thermoplastic elastomer)
Headgear clip, lock	PA (polyamide), POM (polyoxymethylene)
Forehead support, rotating sleeve, mask body, elbow	PA (polyamide)
Headgear	Elastane, polyester, PU (polyurethane), UBL loop, cotton, PA (polyamide)
Anti-asphyxia valve	TPE, PP (polypropylene)
Quick-release cord	POM, polyester, PA (polyamide)

All parts of the mask do not contain latex, PVC (polyvinylchloride) and DEHP (diethylhexyl phthalate).

## 9.2 Characteristic pressure/flow curve

The characteristic pressure/flow curve shows the outlet flow in dependence with therapy pressure.

## 10 Markings and symbols

The following symbols may be applied to the device, the device label, accessories or packaging.

Symbol	Description
	Unique device identifier (uniform product code for medical devices)
	Permitted temperature range for transport and storage
	Use by date
	Keep out of sunlight
	Order number
	Indicates the product is a medical device
	Lot number
	Manufacturer and, if necessary, date of manufacture
	Follow the instructions for use
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

# 11 Warranty

Löwenstein Medical gives the customer a limited manufacturer warranty on new original Löwenstein Medical products and any replacement part fitted by Löwenstein Medical in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the website of the manufacturer. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, contact your specialist dealer.

<b>Product</b>	<b>Warranty periods</b>
Masks including accessories	6 months

# 12 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

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